



CRITIKON CLASSIC-CUF

Blood Pressure Cuffs

CRITIKON™ CLASSIC-CUF™ BP cuffs are designed and clinically validated for limited reuse on multiple patients. Made with durable vinyl construction, this time-tested CRITIKON blood pressure cuff is both easy-to-clean and strong to withstand repeated inflations – yet economical enough for single patient use.



A trusted choice for clinical excellence

- Cuffs are held to the same quality standard as all GE patient monitors, with repeated inflations to ensure accurate readings
- Large bold print helps you quickly choose the right cuff, enabling more accurate measurement

Constructed for patient safety and durability

- Easy-to-clean, vinyl material does not contain natural rubber latex or DEHP, making them safer for your patients
- A full range of cuff sizes are available, from neonate, pediatric, adult and thigh, ensuring that the measurement you receive is from proper fit
- DINACLICK™ connector meets international standards and prevents misconnections to luer fittings

Designed to improve workflow

- Easy-to-use, fast DINACLICK connector enables nurses to link the cuff to the NIBP hose quicker than a pair of traditional tube connectors
- Compliant with AAMI and AHA NIBP cuff sizing standards for prompt, correct sizing
- A full range of adapters and connectors are compatible with automated monitors and a variety of manual devices, including wall manometers and hand-held gauges, enabling a universal connection system throughout the facility

Cuff assortment packs available

- For low-volume cuff utilizers, several types of assortment packs are available for easier ordering and more optimal inventory management



DINACLICK
Connector



CLASSIC-CUF			GE Healthcare				Non-GE Healthcare	
								
Size	Color	Limb Circumference	2-Tube Screw Connector	2-Tube DINACLICK Connector	2-Tube Submin Connector	2-Tube Mated Submin Connector	1-Tube Bayonet Connector	CLASSIC-CUF Inflation System
Infant	Orange	8-13 cm	2618	CLA-P1-2A	2354	2170 ¹	2300	
Child	Green	12-19 cm	2613	CLA-P2-2A	2355	2171	2301	002287
Small Adult	Lt. Blue	17-25 cm	2608	CLA-A1-2A	2356	2172	2302	002285
Adult	Navy	23-33 cm	2603	CLA-A2-2A	2357	2173	2303	002283
Adult Long	Navy	23-33 cm	2647	CLA-A2-2A-L	2352	2174	2306	
Large Adult	Rose	31-40 cm	2643	CLA-A3-2A	2358	2175	2304	002281
Large Adult Long	Rose	31-40 cm	2649	CLA-A3-2A-L	2353	2176	2307	
Thigh	Brown	38-50 cm	2648	CLA-T1-2A	2359	2177	2305	

1. Infant BP Cuff for Datex-Ohmeda* monitors has Male Slip Luer Connector

Note: 20 cuffs/box

CLASSIC-CUF Assortment Pack	Part Number	Description
2-Tube Screw Connectors	2692	1 Infant, 1 Child, 2 Small Adult, 2 Adult, 1 Adult Long, 2 Large Adult, 1 Thigh (10 cuff/box)
2-Tube Submin Connectors	2292	1 Infant, 1 Child, 2 Small Adult, 2 Adult, 1 Adult Long, 2 Large Adult, 1 Thigh (10 cuffs/box)
2-Tube DINACLICK Connectors	2059304-001	5 Small Adult, 10 Adult, 5 Large Adult (20 cuffs/box)
	2059304-002	1 Small Adult, 1 Adult, 1 Large Adult (3 cuffs/pack)

To learn more about GE Healthcare's comprehensive cuff solutions, visit us at cuffs.gehealthcare.com



Imagination at work

www.gehealthcare.com. Product may not be available in all countries and regions. Contact a GE Healthcare Representative for more information. Data subject to change.

© 2017 General Electric Company.

GE, GE Monogram, Imagination at work, CLASSIC-CUF, CRITIKON, and DINALICK are trademarks of General Electric Company.

Reproduction in any form is forbidden without prior written permission from GE. Nothing in this material should be used to diagnose or treat any disease or condition. Readers must consult a healthcare professional.