

Single patient use

Limited reuse

Hospital reusable

General care reusable

YELLOW CRITIKON CLASSIC-CUF

Blood Pressure Cuffs

The yellow CRITIKON™ CLASSIC-CUF™ is designed for single patient use within an isolation environment. These vinyl cuffs are priced as disposable-use cuffs and support your infection prevention protocols. They provide the same quality and assurance that you expect from the entire line of CRITIKON cuffs, while tailored for isolation areas.



A trusted choice for clinical excellence

- Cuffs are held to the same quality standard as all GE patient monitors, with repeated inflations to ensure accurate readings
- Large bold print helps you quickly choose the right cuff, enabling a more accurate measurement
- Yellow color helps make it easy to identify in isolation rooms

Constructed for patient safety and durability


- Easy-to-clean, vinyl material is not made with natural rubber latex or DEHP, making them safer for your patients
- A full range of cuff sizes are available, from infant to adult and thigh, ensuring that the measurement you receive is from a proper fit
- DINACLICK™ connector meets international standards and prevents misconnections to luer fittings

Designed to improve workflow

- Easy-to-use, fast DINACLICK connector enables nurses to link the cuff to the monitor quicker than other dual-tube connectors
- Compliant with AAMI and AHA NIBP cuff sizing standards for prompt, correct sizing
- A full range of adapters and connectors are compatible with automated monitors and a variety of manual devices, including wall manometers and hand-held gauges, enabling a universal connection system throughout the facility



DINACLICK
Connector



YELLOW CLASSIC-CUF

GE Healthcare

Size	Color	Limb Circumference	2-Tube Screw Connector	2-Tube DINACLICK Connector	2-Tube Submin Connector
Infant	Yellow/Blue	8-13 cm	2650	CLA-P1-2A-X	2670
Child	Yellow/Blue	12-19 cm	2651	CLA-P2-2A-X	2671
Small Adult	Yellow/Blue	17-25 cm	2607	CLA-A1-2A-X	2672
Adult	Yellow/Blue	23-33 cm	2602	CLA-A2-2A-X	2673
Large Adult	Yellow/Blue	31-40 cm	2642	CLA-A3-2A-X	2674
Thigh	Yellow/Blue	38-50 cm	2652	CLA-T1-2A-X	2675

Note: 20 cuffs/box.

To learn more about GE Healthcare's comprehensive cuff solutions, visit us at cuffs.gehealthcare.com



Imagination at work

www.gehealthcare.com. Product may not be available in all countries and regions. Contact a GE Healthcare Representative for more information. Data subject to change.

© 2016 General Electric Company.

GE, GE Monogram, Imagination at work, CRITIKON, CLASSIC-CUF, and DINALICK are trademarks of General Electric Company.

Reproduction in any form is forbidden without prior written permission from GE. Nothing in this material should be used to diagnose or treat any disease or condition. Readers must consult a healthcare professional.